


MAE 171A: Heat Transfer Experimental Procedure

Rahul Kapadia


Objectives

- Week 1
 1. To calibrate the heat flux meter
 2. To measure the velocity profile across the duct
- Week 2
 1. Measure the heat flux in forced convection, and free convection
 2. Determine the heat transfer coefficient (h), Nusselt number (Nu). Plot Nusselt numbers (Nu) against Reynolds number (Re) for forced convection, and against Raleigh Number (Ra) for free convection.
- Week 3
 1. Complete any unfinished experiments from the flat plate heat transfer experiment from previous weeks
 2. Determine the heat transfer coefficients for the fins
 3. Plot Nusselt numbers (Nu) against Reynolds Number (Re), or Rayleigh Number (Ra)


Experimental Setup



Control & Measurement of Power



Chanel with Flat Hot Plate



Air flow channel

Flat Hot Plate

Air velocity measurement – Hot
wire anemometer

Finned Hot Plate




Plate temperature
measurement thermocouple

Hot Wire Anemometer



Anemometer Probe

Anemometer Display unit

LabView


TC1 – Measure pre plate /
input Air temperature

TC2 – Plate surface
temperature

TC3 – Measure post
plate / output Air
temperature

Temperature control
knob

Manual / Auto
control knob




Thermocouples

TC3 – Measure post plate /
output Air temperature

TC2 – Plate surface
temperature

TC1 – Measure pre plate /
input Air temperature



Experimental Procedure – Week 1

- Calibration of Heater Power
 - Measure / Record displayed power vs actual power ($P = V * I$)
- Measurement of Velocity Profiles
 - Using Anemometer, measure air velocity at three different locations across the cross-section for 3 different center line velocities (0.6, 1.0, 2.0 m/s)

Experimental Procedure – Week 2: Flat Plate

- Use LabView Auto mode to ramp up plate temperature to within 1-2 °C of desired T, switch to manual and control power to reach desired temperature
- At each temperature, measure actual heater power & air temperatures

Exp #	1	2	3	4	5	6	7	8	9	10	11
Plate Temp (°C)	40	50	60	70	80	60	60	60	60	80	80
Flow speed (m/s)	0	0	0	0	0	0.5	1.0	1.5	2.0	1.0	2.0

Experimental Procedure – Week 3:

Finned Plate

- Change flat plate to Finned plate, and follow same procedure from week 2.
- At each temperature, measure actual heater power & air temperatures

Exp #	1	2	3	4	5	6	7	8	9
Plate Temp (°C)	40	50	60	70	80	40	40	40	40
Flow speed (m/s)	0	0	0	0	0	0.5	1.0	1.5	2.0

Lab Logistics with respect to MAE 170

- No continuous TA / Lab manager supervision
- TA / Lab Manager will discuss at beginning of section, procedure to be performed by students
- When in doubt ask TA / Lab Manager
- If you finish recording data early, try to perform analysis in lab
- Results Analysis more important / time consuming than actual data recording